

DELSONIC 2000 SPARE PARTS LIST

CAT. NO	PART NO.	DESCRIPTION
512101		DELSONIC 2000 NEW HANDPIECE SLEEVE
512102	9,10	DELSONIC 2000 NEW HANDPIECE
512103	9	DELSONIC 2000 NEW SET OF HP SLEEVES (3 PC)
512104	12	DELSONIC 2000 NEW HOLDER COMPLETE
512105	18	DELSONIC 2000 NEW WATER REGULATOR
512106	13,14	DELSONIC 2000 NEW SET OF CONTROL KNOBS
512136	8	DELSONIC 2000 NEW TIP CHANGING TOOL
512001		DELSONIC 2000 PC BOARD
512022	16	DELSONIC 2000 FOOTSWITCH + CABLE
512027		DELSONIC 2000 PACKAGING
506064	1,2	DELSONIC 2000 WATER TUBING 2M+FILTER
147901		DELSONIC 2000 NEW RUBBER FOOT
152101	5	DELSONIC 2000 NEW ON/OFF SWITCH
154004	21	DELSONIC 2000 NEW FUSE 1 AMP.
154007		DELSONIC 2000 NEW PC BOARD 0.05 AMP
158002	17	DELSONIC 2000 NEW CABLE MAINS
312105	19	DELSONIC 2000 WATER SOLENOID
171021	20	DELSONIC 2000 WATER LINE CONNECTOR
171020		

DELSONIC 2000 scaling tips

512050	Tip #1 is a universal tip used for the removal of light to heavy SUPRA-gingival deposits. It can be used on all power settings. It is used for routine scaling procedures.
512052	Tip #2 Round Straight: This tip is designed for removal of light to moderately heavy SUPRA and SUB-gingival deposits. It should not be used on more than medium power setting.
512058	Tip #3 Round Right: . These are indicated for accessing and cleaning open furcation areas. It should not be used on more than medium power setting.
512060	Tip #4 Round Left: . These are indicated for accessing and cleaning open furcation areas. It should not be used on more than medium power setting.
512061	Tip #5 Flat: is also a universal tip used for the removal of light to heavy SUPRA-gingival deposits. It can be used on all power settings. It is used for routine scaling procedures.
512051	Tip # 6 This tip is designed for accessing deep periodontal pockets. It is for removal of light to moderate deposits and should not be used on higher than low to medium power setting.