

**DELDENT MINIBLASTER
FACTS AT A GLANCE
PRODUCT CODE #505000**

- *** FULLY AUTOCLAVABLE CLINICAL SANDBLASTER
PREVENTS CROSS INFECTION**
- *** FEATHERWEIGHT ALUMINUM HEAD & BODY FOR BETTER
CONTROL**
- *** BREAKTHROUGH FINGERTIP-ACTIVATED VALVE AT
FRONT CENTER GIVES MAXIMUM ERGONOMIC
CONTROL FOR BOTH RIGHT AND LEFT-HANDED USERS!**
- *** SPRING-LOADED RACHET SYSTEM SPEEDS
NOZZLE ROTATION AND LOCKS ANGULATION IN A JIFFY**
- *** FDA APPROVED FOR INTRAORAL USE**
- *** BUILT-IN QUICK DISCONNECT ALLOWS SIMPLE REMOVAL
FOR AUTOCLAVING**
- *** USES 50 AND 90 MICRON ALUMINUM OXIDE**
- *** READY TO GO – INCLUDES 8 FT. AIR LINE WITH IN-LINE SHUT
OFF AND PLASTIC T-FITTINGS**

**(ACCESSORIES AVAILABLE, INCLUDING METAL QUICK
DISCONNECT KITS AND EXTENSION NOZZLE TIP FOR HARD-TO-
REACH AREAS WITH INTRAORAL USE)**