

JETSONIC 2000M / Cavitron®JET Plus Comparison chart

JETSONIC 2000M	Cavitron® JET Plus
<b>SCALING TECHNOLOGY</b>	<b>SCALING TECHNOLOGY</b>
Unique Sensor technology For automatic selection of scaling or polishing modes.	
Advanced patented digital circuitry.	
Advance micro-chip feed-back system	
Accepts all 25K and 30K inserts	Only compatible with 30K Inserts
2-position foot control	2-position foot control
Rinse mode available	Rinse mode available
<b>POWER SURGE TECHNOLOGY</b> Double 'clicking' on the foot control activates the Power Surge Mode technology to <b>increase the selected power setting by an additional 15-20%</b> for removal of tenacious deposits. The pre-selected power setting returns automatically when the foot control is released.	<b>SPS TECHNOLOGY</b> Automatically maintains the selected power as the tip comes into contact with the tooth surface
Water heating feature	Water heating feature
<b>POLISHING TECHNOLOGY</b>	<b>POLISHING TECHNOLOGY</b>
<b>H.S.T POLISHING TECHNOLOGY</b> This patented polishing technology generates an homogenous polishing spray that cleans more efficiently	
<b>NON-CLOG POLISHING NOZZEL</b> The single aperture nozzle results in continuous flushing of the nozzle and the prevention of tedious clogging.	
Autoclavable control knobs, handpiece holders and polishing handpiece.	Autoclavable handpiece
<b>KIT CONTENT</b>	<b>KIT CONTENT</b>
1 25K and 1 30K scaling Inserts Polishing handpiece. Scaling handpiece. Replacement pack of autoclavable handpiece holders and control knobs. Medical-grade packaging of polishing powder. Replacement spray head tubing. In-Line water filter Air and water quick- connectors Operating instruction.	1 30K scaling Insert Jet-Mate handpiece Air polishing nozzle 1 bottle of Prophy-Jet powder. 1 bottle of JET-Fresh powder. Foot control and cable 4 pack "AA" batteries. Instruction manual.