

JETSONIC 2000M / Cavitron® JET Plus Comparison Chart

JETSONIC 2000M	Cavitron® Jet Plus
SCALING TECHNOLOGY	SCALING TECHNOLOGY
Unique sensor technology for automatic selection of scaling or polishing modes	
Advanced patented digital circuitry	
Advanced micro-chip feed-back system	
Accepts all 25K and 30K inserts	Only compatible with 30K inserts
2-position foot control	2-position foot control
Rinse Mode available	Rinse Mode available
POWER SURGE TECHNOLOGY Double 'clicking ' on the foot control activates the power surge mode technology to increase the selected power setting by an additional 15-20% for removal of tenacious deposits. The pre- selected power setting returns automatically when the foot control is released	SPS TECHNOLOGY Automatically maintains the selected power as the tip comes into contact with the tooth surface
Water heating feature	Water heating feature
POLISHING TECHNOLOGY	POLISHING TECHNOLOGY
H.S.T. Polishing Technology This patented polishing technology generates an homogenous polishing spray that cleans more efficiently	
NON CLOG POLISHING NOZZEL The single aperture nozzle results in continuous flushing of the nozzle and prevention of tedious clogging	
Autoclavable control knobs, handpiece holders and polishing handpiece	Autoclavable handpiece
KIT CONTENT	KIT CONTENT
1 25K and 1 30K scaling inserts Polishing handpiece Scaling handpiece Replacement pack of autoclavable handpiece holders and control knobs Medical grade packaged polishing powder Replacement spray head tubing In-line water filter Air and Water quick-connectors Operating instructions	1 30K scaling insert Jet –Mate handpiece Air polishing Nozzel 1 bottle of Prophy- Jet powder 1 bottle of JET-Fresh powder Foot control and cable 4 pack "AA" batteries Instruction manual