

Delsonic 2000™ Facts at a glance

- The Delsonic 2000™ combines state of the art piezoceramic technology with a sleek compact design.
- Designed to prevent cross infection it is supplied with 3 autoclavable handpiece sleeves, 2 sets of autoclavable control knobs,
- Advanced circuitry generates linear tip movement at nominal 30 KHz for smoother treated surfaces.
- Automatic tuning and torque adjustment ensures efficient yet gentle calculus removal.
- Scaling tips designed by dentists to meet all clinical needs and to ensure that the water spray reaches the very end of the tip where it is needed even in deep periodontal pockets.
- Supplied with 3 universal scaling tips, 2 tip replacement tools, foot control, water line with in-line water filter and quick disconnect.