



Take a **TIP** from us!



Advanced Dental Technology

DEL DENT

Deldent Ultrasonic inserts give you quality, variety and a price you can't beat !

Power setting, low -moderate.		
25K - DF#3 Flat	508600	
Power setting, low -moderate.		
30K - DF#3 Flat	508700	



- Heavy to moderate supragingival deposits.
- Cement removal.

Power setting, low -moderate.		
25K - DF#10 Universal	508601	
Power setting, low -moderate.		
30K - DF#10 Universal	508701	



- Heavy to moderate supragingival deposits.
- Light to moderate subgingival deposits.

THINLINE SUBGINGIVAL INSERTS		
Power setting, low only.		
25K - DF#100 Straight	508604	
25K - DF#100 Right	508605	
25K - DF#100 Left	508606	
THINLINE SUBGINGIVAL INSERTS		
Power setting, low only.		
30K - DF#100 Straight	508704	
30K - DF#100 Right	508705	
30K - DF#100 Left	508706	



- Light subgingival deposits.

Power setting, low -moderate.		
25K - DF#1000 3-bend	508607	
Power setting, low -moderate.		
30K - DF#1000 3-bend	508707	



- Light to moderate supragingival deposits.

ALL METAL INSERTS - "DP-SERIES"

Power setting, low -moderate.		
25K - DP#10 Universal	508801	
SUBGINGIVAL INSERTS		
Power setting, low only.		
25K - DP#100 Straight	508804	
25K - DP#100 Right	508805	
25K - DP#100 Left	508806	

Power setting, low -moderate.		
30K - DP#10 Universal	508901	
SUBGINGIVAL INSERTS		
Power setting, low only.		
30K - DP#100 Straight	508904	
30K - DP#100 Right	508905	
30K - DP#100 Left	508906	



- DP#10 Heavy to moderate supragingival deposits.
- DP#10 Light to moderate subgingival deposits.
- DP#100 Light subgingival deposits.

DELIDENT ULTRASONIC INSERT COMPARISON CHART *

		Deldent Ltd.		Hu-Friedy ®		Dentsply ®		Whaledent ®
Cat. No.	Handle	Description		Similar to		Similar to		Similar to
508600	Resin	DF#3 Flat	25K	#3 Streamline Beavertail	25K	TFI-3 Thin Beavertail	25K	
508700	Resin	DF#3 Flat	30K	#3 Streamline Beavertail	30K	TFI-3 Thin Beavertail	30K	
508601	Resin	DF#10 Universal	25K	Original Prophy Plus	25K	TFI-10 Straight	25K	#10
508701	Resin	DF#10 Universal	30K	Original Prophy Plus	30K	TFI-10 Straight	30K	#10
508604	Resin	DF#100 Straight	25K	#100 Streamline	25K	Shape similar to FSI-SLI 10	25K	#10 Slim
508605	Resin	DF#100 Right	25K	After Five Right	25K	Shape similar to FSI-SLI 10R	25K	
508606	Resin	DF#100 Left	25K	After Five Left	25K	Shape similar to FSI-SLI 10L	25K	
508704	Resin	DF#100 Straight	30K	#100 Streamline	30K	Shape similar to FSI-SLI 10	30K	#10 Slim
508705	Resin	DF#100 Right	30K	After Five Right	30K	Shape similar to FSI-SLI 10R	30K	
508706	Resin	DF#100 Left	30K	After Five Left	30K	Shape similar to FSI-SLI 10L	30K	
508607	Resin	DF#1000 3-bend	25K	#1000 Streamline Triple-Bend	25K			#1000 Triple-Bend
508707	Resin	DF#1000 3-bend	30K	#1000 Streamline Triple-Bend	30K	TFI-1000 (Triple-Bend)	30K	#1000 Triple-Bend
All Metal Inserts								
508801	Metal	DP#10 Universal	25K	#10 After five	25K	P10	25K	
508901	Metal	DP#10 Universal	30K	#10 After five	30K			
508804	Metal	DP#100 Straight	25K			SLI-10S Straight	25K	
508805	Metal	DP#100 Right	25K			SLI-10R Right	25K	
508806	Metal	DP#100 Left	25K			SLI-10L Left	25K	
508904	Metal	DP#100 Straight	30K			SLI-10S Straight	30K	
508905	Metal	DP#100 Right	30K			SLI-10R Right	30K	
508906	Metal	DP#100 Left	30K			SLI-10L Left	30K	

*This chart is intended as a guide to the user and is based on similarity in shape at the time of printing.

- Finely crafted in Israel
- All-metal and resin handle styles available
- Supreme quality
- The better choice of inserts at highly competitive prices
- Stock up today!

F U L L Y A U T O C L A V A B L E

ULTRASONIC INSERTS
FOR ALL STACK-TYPE
ULTRASONIC SCALER UNITS
SUCH AS
CAVITRON® PARKELL®
AND OTHERS

A selection of Ultrasonic Scaling Inserts which are designed for use on all stack type (magnetostrictive) ultrasonic scaling units such as Cavitron®, Parkell® etc.

‘O’ RING REPLACEMENT PACKS

Cat.No.	Description
500005	Set of 10 Replacement ‘O’ Rings for DF Inserts
500007	Set of 10 Replacement ‘O’ Rings for DP Inserts
500008	Set of 10 Replacement ‘O’ Rings for DF Inserts + Tool
500009	Set of 10 Replacement ‘O’ Rings for DP Inserts + Tool



DELIDENT

Advanced Dental Technology

Deldent Ltd.

19 Keren Kayemet St.

4937226 Petach Tikva, Israel

Tel: 972-3-9049668 Fax: 972-3-9049681

■ e-mail: info@deldent .com ■ www.deldent.com



CE
0344

Your dealer