



**DEL DENT**  
Advanced Dental Technologies

# DiaGlaze™

## Dental Diamond Polishing Paste

Deldent's advanced natural diamond technology has developed a unique polishing system containing ultra-fine, highly graded natural diamond particles that is offered in 3 superbly controlled grades of dental diamond paste.

DiaGlaze, with its pleasant mint flavor, efficiently achieves the ultimate in surface polish for porcelain, gold base metals and composite restorations in minutes.

The 3-part starter kit includes 3 x 2gm. syringes, and optimum results are achieved using a progressive technique, that begins with regular, followed by fine and finished with the extra-fine grade.

Individual refill syringes are available in each of the grades.



C.No. 183090		<b>Deldent Ltd.</b> Petach-Tikva Israel	<b>Regular</b>
		<b>DiaGlaze</b>	
		Diamond Polishing Paste	
		Pasta de Pulido de Diamante	
		Pate a Polir Diamantee	
	CE		
	2gm.	Lot: <input type="text"/>	Fax: 972-3-9049681 WWW: deldent.com

C.No. 183092		<b>Deldent Ltd.</b> Petach-Tikva Israel	<b>Extra Fine</b>
		<b>DiaGlaze</b>	
		Diamond Polishing Paste	
		Pasta de Pulido de Diamante	
		Pate a Polir Diamantee	
	CE		
	2gm.	Lot: <input type="text"/>	Fax: 972-3-9049681 WWW: deldent.com

C.No. 183091		<b>Deldent Ltd.</b> Petach-Tikva Israel	<b>Fine</b>
		<b>DiaGlaze</b>	
		Diamond Polishing Paste	
		Pasta de Pulido de Diamante	
		Pate a Polir Diamantee	
	CE		
	2gm.	Lot: <input type="text"/>	Fax: 972-3-9049681 WWW: deldent.com

### Diaglaze Dental Diamond Paste

- #520004 Starter kit of 3 x 2gm. syringes (Regular, Fine & Extra-fine)
- #520007 Refill - 2gm. Regular grade
- #520006 Refill - 2gm. Fine grade
- #520005 Refill - 2gm. Extra-fine grade

Manufactured by:  
**DEL DENT LTD.**  
19, Keren Kayemet St. Petach Tikva 49372, Israel  
• www.deldent.com • Email: info@deldent.com